

Potassium Sorbate - - Aging Screen

(Conversation w/ Lib Villegas on June 4/14)

Products

6 mg. control (G242 w/ G7-KDN)	PD4684 A
1% Lev. Acid - control #2	PD4684 J
1% Pot. Sorbate (highest level)	PD4684 G

Analyses

Include:

Moisture, Water Activity
Sorbic Acid in tobacco & smoke (6)
Levulinic Acid in " & " (5)

FTC

Sampling / Timing / Conditions

- Ø TIME Thursday, April 16
- Monitor / Sample every 2 weeks for 6 wks.
- Freeze ; 98/20 ; 75/60

Summary

- 3 test parts
- 3 environments
- 4 time points